





Information Questionnaire — Pouch

Customer information	1.Contents			
	Contents	<input type="checkbox"/> Solid	Liquid <input type="checkbox"/> lipidic <input type="checkbox"/> aqueous	<input type="checkbox"/> Powder <input type="checkbox"/> Paste
	Details			
	2.Use Condition			
<input type="checkbox"/> Retort[≥135°C]: ___°C* ___min		<input type="checkbox"/> Semi-retort[121~135°C]: ___°C* ___min		
<input type="checkbox"/> Boiled[80~110°C]: ___°C* ___min		<input type="checkbox"/> Room Temp: ___°C* ___min		
3.Packaging Unit				
<input type="checkbox"/> Manual filling sealed line		<input type="checkbox"/> Auto filling sealed line	<input type="checkbox"/> Others	
4.Products' storage & shelf life				
<input type="checkbox"/> Normal <input type="checkbox"/> Chill <input type="checkbox"/> Freeze <input type="checkbox"/> Other Conditions			Shelf life: ___ days	
Products information	Structure	Foil's brilliant		<input type="checkbox"/> Toward un-sealer <input type="checkbox"/> Toward sealer
	Size[W*L*T]	___mm* ___mm* ___um	Amount	___Pcs
	Printed	<input type="checkbox"/> Yes, Color No.: ___	<input type="checkbox"/> reverse	<input type="checkbox"/> surface <input type="checkbox"/> No
	Packing	<input type="checkbox"/> Cases <input type="checkbox"/> Pallet <input type="checkbox"/> Pallet+case <input type="checkbox"/> Others		
	Pouch Type (Please choose one type and fill the sizes[mm])			
	<input type="checkbox"/> 2 side seal		W: L: S1:	<input type="checkbox"/> 3 side seal
	<input type="checkbox"/> stand up with insert gusset		W: L: S1: S2: G:	<input type="checkbox"/> stand up with folded gusset
	<input type="checkbox"/> stand up (folded gusset/zipper)		W: L: S1: S2: G: Z:	<input type="checkbox"/> 3 side seal with zipper
	<input type="checkbox"/> quad seal		W: L: S1: S2: G:	<input type="checkbox"/> stand up (insert gusset/zipper)
	<input type="checkbox"/> use describe in details below:			<input type="checkbox"/> fin seal
			<input type="checkbox"/> fin seal with side gusset	
			Preparer:	
			Verifier:	